



FOOTWORK néo

THE IDEAL PEDOMETER TO START ITS STATIC, POSTURO AND DYNAMIC ANALYSIS

A TECHNOLOGY SPECIALLY ADAPTED TO PODIATRY

- An indispensable diagnostic tool for your clinical examination.
- The reliability of a real measurement helps you to explore the behavior of the foot
- Quickly and accurately, the functions of the foot are acquired with high quality graphics
- With the Gait Line, the maximum and average pressures, the duration of the supports are essential information to your consultation.
- Compared to the competition, FootWork® is an exceptional value for money with its 2704 capacitive sensors whose accuracy and reliability are of superior quality.

3 NEW DESIGNS TO CHOOSE FROM



NEW DESIGN NEW ELECTRONICS





Inspiration Aluminium

Inspiration Cible

EXAMENS

- Pressure
- Comparison
- Posture
- Tim
- Barefoot / Shoes










STATIC ANALYSIS


Centers of pressure and maximum pressure for each foot, front and rear D/G distribution and by zone, projection of the center of gravity, study zone, film of the static.

FOOTWORK néo®




Stabilometry: ellipses and dimensions of centers of thrust. Oscillation and % in area and time: frontal and sagittal plane. Parameterizable in duration, frequency, % of points.


DYNAMIC ANALYSIS




Average and maximum pressures, pressure duration and pressure/time integral.



Film of the step sequence with simultaneous display of force, surface and pressure curves.



Automatic synchronized video (optional).



Report : Easy and fast editions, sending by mail of the analyses.